

# SmartPeg Reference Guide

Valid from January 2016 | 25016-60

Implant Company	Model	Size	SmartPeg	Art. No
A.B. DENTAL DEVICES	Screw Type, Conical Groovy	3.75	Type 32	100440
	Screw Type, Conical Groovy	4.5	Type 27	100431
ADDENTA	Bifasico	4.0	Type 14	100374
ADIN DENTAL IMPLANT SYSTEM	One	3.3	Type 17	100390
	One	4.2	Type 17	100390
	Swell	3.3	Type 49	100480
	Swell	3.75	Type 49	100480
	Swell	4.2	Type 49	100480
	Swell	5.0	Type 49	100480
	Swell	6.0	Type 49	100480
	Touareg CloseFit	NP	Type 62	100545
	Touareg CloseFit	RP	Type 60	100543
	Touareg CloseFit	WP	Type 49	100480
	Touareg-S	4.2	Type 49	100480
	Touareg-S	5.0	Type 49	100480
ALPHA BIO IMPLANT	ATI		Type 26	100425
	ATID	3.75	Type 26	100425
	SPI		Type 26	100425
ALPHA-GATE	Bioactive	3.75	Type 35	100446
ALPHATECH		3.4	Type 46	100477
		3.8	Type 46	100477
		4.3	Type 46	100477
		5.0	Type 46	100477
		6.0	Type 46	100477
	7.0	Type 46	100477	
ANKYLOS	C/X	A	Type 16	100388
	C/X	B	Type 16	100388
	C/X	C	Type 16	100388
ANTHOGYR	Anthofit External Hexagon HE		Type 01	100353
	Anthofit Internal Octagon OI		Type 18	100398
	Axiom PX	3.4	Type 38	100455
	Axiom PX	4.0	Type 38	100455
	Axiom PX	4.6	Type 38	100455
	Axiom PX	5.2	Type 38	100455
	Axiom REG	3.4	Type 38	100455
Axiom REG	4.0	Type 38	100455	

Implant Company	Model	Size	SmartPeg	Art. No
<b>ANTHOGYR (Continued)</b>	Axiom REG	4.6	Type 38	100455
	Axiom REG	5.2	Type 38	100455
	Ossfit		Type 04	100350
<b>ARDS IMPLANTS</b>	Classic		Type 27	100431
	Smart		Type 27	100431
	Premium		Type 27	100431
<b>ASTRA TECH (Please see DENTSPLY)</b>				
<b>AVINENT IMPLANT SYSTEM</b>	Coral External	3.5	Type 30	100436
	Coral External	4.1	Type 01	100353
	Coral External	5.1	Type 03	100357
	Ocean External	3.5	Type 30	100436
	Ocean External	4.1	Type 01	100353
	Internal Hexagon Coral	3.5	Type 32	100440
	Internal Hexagon Coral	4.1	Type 32	100440
	Internal Hexagon Ocean	3.5	Type 32	100440
	Internal Hexagon Ocean	4.1	Type 32	100440
<b>AXIS</b>	Ceramic	4.0	Type 38	100455
<b>BEGO IMPLANT SYSTEMS</b>	Semados	RI 3.75	Type 26	100425
	Semados	RI 4.1	Type 26	100425
	Semados	RI 4.5	Type 26	100425
	Semados	RI 5.5	Type 26	100425
	Semados	S 3.25	Type 26	100425
	Semados	S 3.75	Type 26	100425
	Semados	S 4.1	Type 26	100425
	Semados	S 4.5	Type 26	100425
	Semados	S 5.5	Type 26	100425
	Semados Mini	2.7	Type 02	100355
	Semados Mini	2.9	Type 02	100355
	Semados Mini	3.1	Type 02	100355
	Semados RS/R SX	3.0	Type 50	100481
	Semados RS/R SX	3.75	Type 26	100425
	Semados RS/R SX	4.1	Type 26	100425
	Semados RS/R SX	4.5	Type 26	100425
	Semados RS/R SX	5.5	Type 26	100425
<b>BIODENTA</b>	Tissue level	NP	Type 21	100402
	Tissue level	RP	Type 47	100478
	Tissue level	WP	Type 47	100478
<b>BIOHORIZONS</b>	External	3.5	Type 30	100436
	External	4.0	Type 01	100353
	External	5.0	Type 01	100353

Implant Company	Model	Size	SmartPeg	Art. No
<b>BIOHORIZONS (Continued)</b>	External	6.0	Type 01	100353
	Internal	Blue platform	Type 33	100442
	Internal	Green platform	Type 27	100431
	Internal	Yellow platform	Type 32	100440
	Internal Laser Lok	3.0	Type 26	100425
	Single Stage	Blue platform	Type 33	100442
	Single Stage	Green platform	Type 27	100431
	Single Stage	Yellow platform	Type 32	100440
	Tapered Internal	Blue platform	Type 33	100442
	Tapered Internal	Green platform	Type 27	100431
	Tapered Internal	Yellow platform	Type 32	100440

<b>BIOMET 3i</b>	Osseotite Certain	XP 4/5	Type 15	100386
	Osseotite Certain	XP 5/6	Type 15	100386
	Osseotite Certain Microminiplant	3.25	Type 45	100472
	Nanotite Certain Microminiplant	3.25	Type 45	100472
	Osseotite External Hex	3.75	Type 01	100353
	Osseotite External Hex	4.0	Type 01	100353
	Osseotite External Hex	5.0	Type 01	100353
	Osseotite External Hex	6.0	Type 01	100353
	Osseotite External Hex	XP 3/4	Type 01	100353
	Osseotite External Hex	XP 4/5	Type 01	100353
	Osseotite External Hex	XP 5/6	Type 01	100353
	Nanotite External Hex	3.75	Type 01	100353
	Nanotite External Hex	4.0	Type 01	100353
	Nanotite External Hex	5.0	Type 01	100353
	Nanotite External Hex	6.0	Type 01	100353
	Nanotite External Hex	XP 3/4	Type 01	100353
	Nanotite External Hex	XP 4/5	Type 01	100353
	Nanotite External Hex	XP 5/6	Type 01	100353
	Osseotite External Hex MicroMiniplant	3.25	Type 01	100353
	Nanotite External Hex MicroMiniplant	3.25	Type 01	100353
	Osseotite Tapered Certain	3.25	Type 45	100472
	Osseotite Tapered Certain	4.0	Type 15	100386
	Osseotite Tapered Certain	5.0	Type 15	100386
	Osseotite Tapered Certain	6.0	Type 15	100386
	Nanotite Tapered Certain	3.25	Type 45	100472
	Nanotite Tapered Certain	4.0	Type 15	100386
	Nanotite Tapered Certain	5.0	Type 15	100386
	Nanotite Tapered Certain	6.0	Type 15	100386
	Osseotite Tapered Certain Prevail	4x3	Type 45	100472
	Osseotite Tapered Certain Prevail	4x5	Type 15	100386
	Osseotite Tapered Certain Prevail	5x6	Type 15	100386
	Nanotite Tapered Certain Prevail	4x3	Type 45	100472
	Nanotite Tapered Certain Prevail	4x5	Type 15	100386
	Nanotite Tapered Certain Prevail	5x6	Type 15	100386
	T3 Implant Internal Connection		Type 15	100386

Implant Company		Model	Size	SmartPeg	Art. No
<b>BIOMET 3i</b> (Continued)		T3 Implant External Connection		Type 01	100353
<b>BIOTEC SRL</b>		EN	All	Type 30	100436
		ER	All	Type 01	100353
		KR	3.05 platf	Type 38	100455
		KW	3.65 platf	Type 32	100440
		IM	All	Type 27	100431
		IR	All	Type 32	100440
		SR	All	Type 06	100378
		TWO	All	Type 32	100440
<b>BONE SYSTEM</b>		2P	5.0	Type 36	100447
<b>BRAINBASE</b>		Brainbase	C	Type 47	100478
		Mytis Arrow Implants	B	Type 04	100350
		Mytis Arrow Implants	BW	Type 04	100350
		Mytis Arrow Implants	E	Type 04	100350
		Mytis Arrow Implants	EW	Type 04	100350
<b>BREIDENT</b>		blueSKY		Type 49	100480
<b>BTI BIOTECHNOLOGY INSTITUTE</b>		Externa	Universal	Type 01	100353
		Externa	Universal plus	Type 01	100353
		Interna	Universal	Type 27	100431
		Interna	Universal plus	Type 27	100431
		Interna	Wide	Type 27	100431
		Tiny	2.5	Type 09	100368
		Tiny	3.75	Type 09	100368
<b>BTLOCK</b>		CV1	3.75	Type 49	100480
		CV1	4.5	Type 49	100480
		CV1	5.5	Type 49	100480
		CV1	6.5	Type 49	100480
		CV2	3.75	Type 49	100480
		CV2	4.5	Type 49	100480
		CV2	5.5	Type 49	100480
		CV2	6.5	Type 49	100480
		CV3	3.75	Type 49	100480
		CV3	4.5	Type 49	100480
		CV3	5.5	Type 49	100480
		CV3	6.5	Type 49	100480
<b>CAMLOG</b>		Camlog Cylinder-line	3.3	Type 28	100432
		Camlog Cylinder-line	3.8	Type 23	100411
		Camlog Cylinder-line	4.3	Type 23	100411
		Camlog Cylinder-line	5.0	Type 24	100412
		Camlog Cylinder-line	6.0	Type 24	100412

Implant Company	Model	Size	SmartPeg	Art. No
<b>CAMLOG (Continued)</b>	Camlog Root-line Promote	3.8	Type 23	100411
	Camlog Root-line Promote	4.3	Type 23	100411
	Camlog Root-line Promote	5.0	Type 24	100412
	Camlog Root-line Promote	6.0	Type 24	100412
	Camlog Screw-line	3.3	Type 28	100432
	Camlog Screw-line	3.8	Type 23	100411
	Camlog Screw-line	4.3	Type 23	100411
	Camlog Screw-line	5.0	Type 24	100412
	Camlog Screw-line	6.0	Type 24	100412
	Conelog Screw-line	3.3	Type 65	100557
	Conelog Screw-line	3.8	Type 58	100541
	Conelog Screw-line	4.3	Type 58	100541
	Conelog Screw-line	5.0	Type 07	100380
	iSy (direct to implant)	All	Type 57	100497
	iSy (to implant base)	All	Type 35	100446
<b>CHAORUM</b>	External	3.3	Type 01	100353
	Internal	4.0	Type 47	100478
	Internal	4.8	Type 47	100478
<b>COCHLEAR</b>	Baha BI 300	3.0	Type 30	100436
	Baha BI 300	4.0	Type 30	100436
	Baha BIA 300	6.0	Type 55	100487
	Baha BIA 300	9.0	Type 55	100487
	Baha BIA 400	6.0	Type 55	100487
	Baha BIA 400	8.0	Type 55	100487
	Baha BIA 400	10.0	Type 55	100487
	Baha BIA 400	12.0	Type 55	100487
	Baha BIA 400	14.0	Type 55	100487
	Baha & Vistafix Flange Fixture Hexagonal Connection	3.0	Type 09	100368
	Baha & Vistafix Flange Fixture Hexagonal Connection	4.0	Type 09	100368
	Vistafix VXI 300	3.0	Type 30	100436
	Vistafix VXI 300	4.0	Type 30	100436
<b>COMPRESS</b>		4.0	Type 02	100355
<b>CORTEX</b>	Classix	3.3	Type 26	100425
	Classix	3.8	Type 26	100425
	Classix	4.2	Type 26	100425
	Classix	5.0	Type 26	100425
	Classix	6.0	Type 26	100425
	Dynamix	3.3	Type 26	100425
	Dynamix	3.8	Type 26	100425
	Dynamix	4.2	Type 26	100425
	Dynamix	5.0	Type 26	100425
	Dynamix	6.0	Type 26	100425
	Saturn	3.3	Type 26	100425

Implant Company	Model	Size	SmartPeg	Art. No
<b>CORTEX (Continued)</b>	Saturn	3.8	Type 26	100425
	Saturn	4.2	Type 26	100425
	Saturn	5.0	Type 26	100425
	Saturn	6.0	Type 26	100425
<b>COWELLMEDI</b>	Atlas External	3.3	Type 02	100355
	Atlas External	4.1	Type 01	100353
	Atlas External	5.0	Type 01	100353
	Atlas Submerged	4.0	Type 21	100402
<b>C-TECH</b>	EL	3.5	Type 26	100425
	EL	4.3	Type 26	100425
	EL	5.1	Type 26	100425
<b>DENTAL TECH</b>		WP	Type 47	100478
<b>DENTAURUM</b>	Tiologic	S	Type 45	100472
	Tiologic	M	Type 45	100472
	Tiologic	L	Type 45	100472
<b>DENTEGRIS</b>	Straight & Tapered	3.3	Type 26	100425
	Straight & Tapered	3.75	Type 32	100440
	Straight & Tapered	4.5	Type 32	100440
	Straight & Tapered	4.75	Type 32	100440
	Straight & Tapered	5.5	Type 47	100478
<b>DENTIS DENTAL</b>	Internal	3.7	Type 04	100350
	Internal	4.3	Type 04	100350
	Internal	4.8	Type 04	100350
	OneQ S-clean Narrow (SLA surface/Submerged)	3.0	Type 22	100404
	OneQ S-clean Narrow (SLA surface/Submerged)	3.3	Type 22	100404
	OneQ S-clean Regular (SLA surface/Submerged)	3.7	Type 11	100372
	OneQ S-clean Regular (SLA surface/Submerged)	4.2	Type 11	100372
	OneQ S-clean Regular (SLA surface/Submerged)	4.7	Type 11	100372
	OneQ S-clean Regular (SLA surface/Submerged)	5.2	Type 11	100372
	OneQ S-clean Wide (SLA surface/Submerged)	6.0	Type 11	100372
	OneQ S-clean Wide (SLA surface/Submerged)	7.0	Type 11	100372
	OneQ S-clean Wide (SLA surface/Submerged)	8.0	Type 11	100372
	Save Submerged	5.5	Type 07	100380
	Save Submerged	6.0	Type 07	100380
	S-clean Mini (Submerged)	3.7	Type 11	100372
	S-clean Regular (Submerged)	4.1	Type 11	100372
	S-clean Regular (Submerged)	4.3	Type 11	100372
	S-clean Wide (Submerged)	4.8	Type 11	100372

Implant Company	Model	Size	SmartPeg	Art. No
<b>DENTIUM</b>	Genoss	3.0	Type 63	100554
	Genoss	3.6	Type 21	100402
	Genoss	4.3	Type 21	100402
	Genoss	5.0	Type 21	100402
	Implantium	3.4 (3.6)	Type 06	100378
	Implantium	3.8 (4.0)	Type 06	100378
	Implantium	4.3 (4.5)	Type 06	100378
	Simple Line	3.8	Type 07	100380
	Simple Line	4.3	Type 07	100380
	Simple Line	4.8	Type 07	100380
	Super Line	4.0	Type 07	100380
	Super Line	7.0	Type 07	100380

<b>DENTSPLY IMPLANTS</b>	ASTRA TECH Implant system MicroThread	3.5	Type 05	100376	
	ASTRA TECH Implant system MicroThread	4.0	Type 06	100378	
	ASTRA TECH Implant system MicroThread	4.5	Type 07	100380	
	ASTRA TECH Implant system MicroThread	5.0	Type 07	100380	
	ASTRA TECH Implant system OsseoSpeed Yellow	3.0S	Type 50	100481	
	ASTRA TECH Implant system OsseoSpeed Aqua	3.5S	Type 38	100455	
	ASTRA TECH Implant system OsseoSpeed Aqua	4.0S	Type 38	100455	
	ASTRA TECH Implant system OsseoSpeed Lilac	4.5	Type 07	100380	
	ASTRA TECH Implant system OsseoSpeed Lilac	5.0	Type 07	100380	
	ASTRA TECH Implant system OsseoSpeed Lilac	5.0S	Type 07	100380	
	ASTRA TECH Implant system OsseoSpeed TX Yellow	3.0S	Type 50	100481	
	ASTRA TECH Implant system OsseoSpeed TX Aqua	3.5S	Type 38	100455	
	ASTRA TECH Implant system OsseoSpeed TX Aqua	4.0S	Type 38	100455	
	ASTRA TECH Implant system OsseoSpeed TX Lilac	4.5	Type 07	100380	
	ASTRA TECH Implant system OsseoSpeed TX Lilac	5.0	Type 07	100380	
	ASTRA TECH Implant system OsseoSpeed TX Lilac	5.0S	Type 07	100380	
	ASTRA TECH Implant system OsseoSpeed TX Profile	4.5	Type 06	100378	
	ASTRA TECH Implant system OsseoSpeed TX Profile	5.0	Type 06	100378	
	ASTRA TECH Implant system OsseoSpeed TX Profile	5.0S	Type 06	100378	
	ASTRA TECH Implant system OsseoSpeed EV	3.0S Green	Type 50	100481	
	ASTRA TECH Implant system OsseoSpeed EV	3.6S Purple	Type 38	100455	
	ASTRA TECH Implant system OsseoSpeed EV	4.2S Yellow	Type 49	100480	
	ASTRA TECH Implant system OsseoSpeed EV	4.8S Blue	Type 47	100478	
	ASTRA TECH Implant system OsseoSpeed EV	5.4S Brown	Type 40	100462	
	ASTRA TECH Implant system OsseoSpeed EV	4.2 PC	Type 10	100370	
	ASTRA TECH Implant system OsseoSpeed EV	4.8 PC	Type 32	100440	
	ASTRA TECH Implant system OsseoSpeed EV	4.2 PS	Type 10	100370	
	ASTRA TECH Implant system OsseoSpeed EV	4.8 PS	Type 32	100440	
		ANKYLOS C/X Implants	A	Type 16	100388
		ANKYLOS C/X Implants	B	Type 16	100388
		ANKYLOS C/X Implants	C	Type 16	100388
		XiVe	3.0	Type 22	100404
		XiVe	3.4	Type 21	100402
		XiVe	3.8	Type 45	100472

Implant Company	Model	Size	SmartPeg	Art. No
<b>DENTSPLY IMPLANTS</b> (Continued)		XiVe 4.5	Type 45	100472
		XiVe 5.5	Type 45	100472
		XiVe TG	Type 08	100366
<b>DIO</b>		External FSN 3.3	Type 02	100355
		External FSN 3.75	Type 01	100353
		External FTN 5.0	Type 03	100357
		SM Submerged 3.8	Type 49	100480
		SM Submerged 4.5	Type 47	100478
		SM Submerged 5.3	Type 47	100478
<b>DTI</b>		3.5	Type 57	100497
		4.2	Type 07	100380
<b>DYNA DENTAL</b>		Dyna Helix DC 3.2	Type 27	100431
		Dyna Helix DC 3.6	Type 27	100431
		Dyna Helix DC 4.2	Type 27	100431
		Dyna Helix DC 5.0	Type 27	100431
		Dyna Helix ST 3.6	Type 27	100431
		Dyna Helix ST 4.2	Type 27	100431
		Dyna Helix ST 5.0	Type 27	100431
<b>ECKERMANN</b>		All Spiral	Type 17	100390
		Hexagono Ext Hex	Type 01	100353
		Triplo	Type 49	100480
<b>EMFILS</b>		Colosso IH ECI Evolution 3.3 x 16	Type 11	100372
		Colosso IH ECI Evolution 4.0 x 7	Type 11	100372
		Colosso IH ECI Evolution 5.0 x 16	Type 11	100372
<b>ENTIFIC</b>		Flange Fixture	Type 01	100353
<b>EUROTEKNIKA</b>		Aesthetica+2 3.6/4.2	Type 61	100544
		Aesthetica+2 3.6/4.8	Type 61	100544
		Aesthetica+2 4.1/4.2	Type 61	100544
		Aesthetica+2 4.1/4.8	Type 61	100544
		Aesthetica+2 4.8/4.8	Type 61	100544
		Aesthetica+2 4.8/6.5	Type 61	100544
		Naturall 3.5	Type 38	100455
		Naturall 4.0	Type 06	100378
		Naturall 4.5	Type 06	100378
		Naturall 5.0	Type 06	100378
		Uneva+ 3.6	Type 01	100353
		Uneva+ 4.1	Type 01	100353
		Uneva+ 4.8	Type 01	100353
		Uneva+ 6.0	Type 01	100353



Implant Company	Model	Size	SmartPeg	Art. No
<b>FAIR IMPLANT</b>	FairTwo	3.5	Type 53	100484
	FairTwo	4.2	Type 53	100484
	FairTwo	5.0	Type 06	100378
	FairTwo	6.0	Type 06	100378
<b>GC IMPLANT</b>	Aadva	Narrow	Type 38	100455
	Aadva	Regular	Type 58	100541
	Aadva	Wide	Type 58	100541
	External Plus	3.0	Type 21	100402
	External Plus	3.4	Type 21	100402
	External Plus	3.8	Type 01	100353
	External Plus	5.0	Type 01	100353
	GENESiO		Type 21	100402
	Internal Plus	3.0	Type 21	100402
	Internal Plus	3.4	Type 21	100402
	SETiO		Type 01	100353
<b>GHIMAS</b>	PFON EVO SERIES	3.6	Type 38	100455
	PFON EVO SERIES	4.0	Type 38	100455
	PFON EVO SERIES	5.0	Type 07	100380
<b>HIOSSEN</b>	ETIII	4.0	Type 47	100478
	ETIII	4.5	Type 47	100478
	ETIII	5.0	Type 47	100478
	ETIII	Mini 3.5	Type 05	100376
	ETIII	UW 6.0	Type 06	100378
	ETIII	UW 7.0	Type 06	100378
	HGII	Mini	Type 05	100376
	HGII	Regular	Type 06	100378
	HGII	Ultra Wide	Type 06	100378
	HGII	Wide	Type 06	100378
	HGIII	Mini	Type 05	100376
HSII	Regular	Type 04	100350	
HSII	Wide	Type 04	100350	
<b>IBS INNO BIO SURG</b>	Magicore	5.5	Type 38	100455
<b>IDI EVOLUTION</b>	Vite	3.75	Type 18	100398
	Vite	4.0	Type 18	100398
<b>IMPLANT DIRECT</b>	Legacy3	3.2	Type 38	100455
	ScrewPlant	3.7	Type 26	100425
	ScrewPlant	4.7	Type 27	100431
	ScrewPlant	5.7	Type 27	100431
<b>IMTEC</b>	Hexed-head	3.75	Type 01	100353
	Hexed-head	4.75	Type 01	100353

Implant Company	Model	Size	SmartPeg	Art. No
<b>INTRA-LOCK SYSTEM</b>	Internal	SQ	Type 06	100378
	Internal	Standard	Type 47	100478
	Internal	Wide	Type 47	100478
<b>ITERUM</b>	Conical	3.75	Type 26	100425
<b>JDENTALCARE</b>	JDEvolution	3.2	Type 38	100455
	JDEvolution	3.7	Type 32	100440
	JDEvolution	4.3	Type 32	100440
	JDEvolution	5.0	Type 32	100440
<b>JMM</b>	POI 2 piece	3.2	Type 18	100398
	POI 2 piece	3.7	Type 18	100398
	POI 2 piece	4.2	Type 18	100398
	POIEX	3.7	Type 15	100386
	POIEX	4.2	Type 49	100480
	POIEX	4.7	Type 47	100478
	POIEX	5.2	Type 47	100478
<b>KENTEC</b>	SB1	3.3	Type 06	100378
	SB2	4.2	Type 06	100378
<b>KEYSTONE DENTAL</b>	Genesis	3.8	Type 27	100431
	Genesis	4.5	Type 33	100442
	Genesis	5.5	Type 33	100442
	Genesis	6.5	Type 33	100442
	Prima Connex	RD 4.1	Type 33	100442
	Prima Connex	SD 3.5	Type 27	100431
	Prima Connex	WD 5.0	Type 33	100442
	Renova RDI	4.5	Type 27	100431
	Renova SDI	3.75	Type 26	100425
	Restore	3.3	Type 30	100436
	Restore	3.75	Type 01	100353
	Restore	4.0	Type 01	100353
	Restore	5.0	Type 03	100357
	Restore	6.0	Type 03	100357
	<b>KLOCKNER IMPLANT SYSTEM</b>	Essential Cone	3.5	Type 49
Essential Cone		4.0	Type 49	100480
Essential Cone		4.5	Type 49	100480
Essential Cone		4.8	Type 49	100480
NK2		3.8	Type 30	100436
S3M			Type 30	100436
S4			Type 30	100436
S6			Type 30	100436
SK			Type 30	100436
SK2			Type 30	100436

Implant Company	Model	Size	SmartPeg	Art. No
<b>KLOCKNER IMPLANT SYSTEM</b>	KL Implant NP	3.5	Type 02	100355
(Continued)	KL Implant NP	3.7	Type 01	100353
	KL Implant NP	4.2	Type 01	100353
	KL Implant WP	5.1	Type 03	100357
	MV VEGA Implant	3.0	Type 57	100497
	NV VEGA Implant	3.5	Type 57	100497
	RV VEGA Implant	4.0	Type 26	100425
	RV VEGA Implant	4.5	Type 26	100425
<b>LASAK</b>	BioniQ Q-lock	Regular	Type 57	100497
	Impladent	2.9	Type 02	100355
	Impladent	3.7	Type 36	100447
	Impladent	5.0	Type 37	100448
<b>LEADER</b>	Implus EH	3.3	Type 01	100353
	Implus EH	3.75	Type 01	100353
	Implus EH	5.0	Type 01	100353
	Implus IH	3.3	Type 26	100425
	Tixos IH	3.75	Type 26	100425
	Tixos IH	4.5	Type 26	100425
	Tixos IH	5.5	Type 26	100425
	Tixos ITMC	3.75	Type 57	100497
	Tixos ITMC	4.5	Type 57	100497
	Tixos ITMC	5.5	Type 57	100497
<b>MEDENTIKA</b>	M-Implant®	3.5	Type 38	100455
	M-Implant®	4.0	Type 38	100455
	M-Implant®	4.5	Type 38	100455
	M-Implant®	5.0	Type 38	100455
<b>MEDICAL INSTINCT</b>	Bone Trust	3.4	Type 44	100466
	Bone Trust	4.0	Type 44	100466
	Bone Trust	5.0	Type 44	100466
	Bone Trust Plus Cone	3.4	Type 44	100466
	Bone Trust Plus Cone	4.0	Type 44	100466
<b>MEDIGMA TECHNOLOGIES</b>	Fix 1		Type 17	100390
	Fix-A-Dent		Type 49	100480
	FixTite		Type 49	100480
	FixTite-S		Type 49	100480
<b>MEGAGEN</b>	AnyOne	3.5	Type 07	100380
	AnyOne	4.0	Type 07	100380
	AnyOne External	4.5	Type 04	100350
	Anyridge		Type 27	100431
	EZ Plus Internal	Small	Type 38	100455
	EZ Plus Internal	Regular	Type 07	100380

Implant Company	Model	Size	SmartPeg	Art. No
<b>MEGAGEN</b> (Continued)	EZ Plus Internal	Wide	Type 07	100380
	EZ Plus Internal	4.0	Type 07	100380
	EZ Plus Internal	4.5	Type 07	100380
	EZ Plus Internal	5.0	Type 07	100380
	Mini Internal	3.0	Type 50	100481
	Resque External	7.0	Type 03	100357
<b>MICRODENT</b>	Duraplant 2.2	4.1	Type 10	100370
	Microdent system		Type 01	100353
	MK	5.6	Type 03	100357
	MS	3.5	Type 30	100436
	MST	3.5	Type 30	100436
	Trylogic	3.8	Type 27	100431
	Trylogic	4.3	Type 27	100431
	Universal	4.1	Type 01	100353
	Universal	5.1	Type 01	100353
<b>MIS</b>	C1 (conical connection)	Standard 3.75	Type 49	100480
	Internal Hexagon (seven, biocom)	3.75	Type 32	100440
	Internal Octagon (mistral)	4.75	Type 04	100350
	Seven Narrow	3.3	Type 63	100554
<b>MOZO-GRAU</b>	MG-Inhex	Std	Type 38	100455
	MG-Inhex	Maxi	Type 07	100380
	MG-Osseous	Mini	Type 01	100353
	MG-Osseous	Std	Type 01	100353
	MG-Osseous	Maxi	Type 03	100357
<b>NEMRIS</b>	Aesthura Classic	3.75	Type 35	100446
	Aesthura Classic	4.75	Type 35	100446
<b>NEOBIOTECH</b>	Neoplant External	4.0	Type 01	100353
	Neoplant External	5.0	Type 01	100353
<b>NEODENT</b>	Alvim CM	3.5	Type 16	100388
	Alvim CM	4.3	Type 16	100388
	Alvim CM	5.0	Type 16	100388
	Alvim II Plus	3.75	Type 21	100402
	Alvim II Plus	4.3	Type 21	100402
	Alvim II Plus	5.0	Type 21	100402
	Alvim Ti	4.3	Type 01	100353
	Titamax CM	3.5	Type 16	100388
	Titamax CM	3.75	Type 16	100388
	Titamax CM	4.0	Type 16	100388
	Titamax CM	5.0	Type 16	100388
	Titamax GT	3.7	Type 09	100368
	Titamax GT	4.0	Type 09	100368

Implant Company	Model	Size	SmartPeg	Art. No
<b>NEODENT</b> (Continued)	Titamax GT	4.5	Type 09	100368
	Titamax II Plus	3.75	Type 21	100402
	Titamax II Plus	5.0	Type 21	100402
	Titamax Ti	3.75	Type 01	100353
	Titamax Ti	4.0	Type 01	100353
	Titamax Ti	4.5	Type 03	100357
	Titamax Ti	5.0	Type 03	100357
<b>NEOSS</b>	Neoss Implant	3.25	Type 38	100455
	ProActive Straight	3.5	Type 18	100398
	ProActive Straight	4.0	Type 18	100398
	ProActive Straight	4.5	Type 18	100398
	ProActive Straight	5.0	Type 18	100398
	ProActive Straight	5.5	Type 18	100398
	ProActive Straight	6.0	Type 18	100398
	ProActive Tapered	3.5	Type 18	100398
	ProActive Tapered	4.0	Type 18	100398
	ProActive Tapered	4.5	Type 18	100398
	ProActive Tapered	5.0	Type 18	100398
	ProActive Tapered	5.5	Type 18	100398
	ProActive Tapered	6.0	Type 18	100398
<b>NOBEL BIOCARE</b>	Brånemark system	NP	Type 02	100355
	Brånemark system	RP	Type 01	100353
	Brånemark system	WP	Type 03	100357
	Brånemark system Novum		Type 01	100353
	Brånemark system Shorty	RP	Type 01	100353
	Brånemark system Shorty	WP	Type 03	100357
	NobelActive	NP	Type 60	100543
	NobelActive	RP	Type 61	100544
	NobelActive	WP	Type 61	100544
	NobelActive	WP 5.5	Type 04	100350
	NobelActive	3.0	Type 62	100545
	NobelDirect	NP	Type 08	100366
	NobelDirect	RP	Type 09	100368
	NobelDirect	WP	Type 09	100368
	NobelParallel Conical Connection	NP	Type 60	100543
	NobelParallel Conical Connection	RP	Type 61	100544
	NobelParallel Conical Connection	WP	Type 61	100544
	NobelReplace Conical Connection	NP	Type 60	100543
	NobelReplace Conical Connection	RP	Type 61	100544
	NobelReplace Conical Connection PMC	NP	Type 60	100543
	NobelReplace Conical Connection PMC	RP	Type 61	100544
	NobelReplace Platform Shift	NP	Type 12	100359
	NobelReplace Platform Shift	RP	Type 13	100361
	NobelReplace Platform Shift	WP	Type 13	100361
	NobelReplace Straight Groovy	NP	Type 12	100359

Implant Company	Model	Size	SmartPeg	Art. No
<b>NOBEL BIOCARE (Continued)</b>	NobelReplace Straight Groovy	RP	Type 13	100361
	NobelReplace Straight Groovy	WP	Type 13	100361
	NobelReplace Tapered Groovy	NP	Type 12	100359
	NobelReplace Tapered Groovy	RP	Type 13	100361
	NobelReplace Tapered Groovy	WP	Type 13	100361
	NobelReplace Tapered Groovy	6.0	Type 40	100462
	NobelSpeedy Groovy	NP	Type 02	100355
	NobelSpeedy Groovy	RP	Type 01	100353
	NobelSpeedy Groovy	WP	Type 03	100357
	NobelSpeedy Replace	NP	Type 12	100359
	NobelSpeedy Replace	RP	Type 13	100361
	NobelSpeedy Replace	WP	Type 13	100361
	NobelSpeedy Replace	6.0	Type 40	100462
	NobelSpeedy Shorty	RP	Type 01	100353
	NobelSpeedy Shorty	WP	Type 03	100357
	Replace Select Straight	NP	Type 12	100359
	Replace Select Straight	RP	Type 13	100361
	Replace Select Straight	WP	Type 13	100361
	Replace Select Tapered	NP	Type 12	100359
	Replace Select Tapered	RP	Type 13	100361
	Replace Select Tapered	WP	Type 13	100361
	Replace Select Tapered	6.0	Type 40	100462
	Replace Select Tapered PMC	NP	Type 12	100359
	Replace Select Tapered PMC	RP	Type 13	100361
	Replace Select Tapered PMC	WP	Type 13	100361
	Replace Select Tapered PMC	6.0	Type 40	100462
	Replace Select TC	NP	Type 12	100359
	Replace Select TC	RP	Type 13	100361
<b>NUCLEOSS</b>	Nucleoss T3	3.4	Type 05	100376
	Nucleoss T3	3.8	Type 05	100376
	Nucleoss T3	4.2	Type 05	100376
	Nucleoss T3	5.0	Type 05	100376
	Nucleoss T4	3.4	Type 05	100376
	Nucleoss T4	3.8	Type 05	100376
	Nucleoss T4	4.2	Type 05	100376
	Nucleoss T5	5.0	Type 07	100380
<b>OCO BIOMEDICAL</b>	Enagage	3.2	Type 38	100455
	Enagage	3.7/3.5	Type 26	100425
	Enagage	4.7/4.5	Type 26	100425
	Enagage	5.7/4.5	Type 27	100431
	ERI	4.0	Type 07	100380
	ERI	5.0	Type 07	100380
	Macro	6.0	Type 47	100478
	TSI	3.25	Type 35	100446
	TSI	4.0	Type 07	100380

Implant Company	Model	Size	SmartPeg	Art. No
<b>OCO BIOMEDICAL</b> (Continued)	TSI	5.0	Type 07	100380
<b>OSSEOLINK</b>		3.5	Type 21	100402
		4.0	Type 21	100402
		5.0	Type 21	100402
<b>OSSTEM</b>	GSII	Mini	Type 05	100376
	GSII	Regular	Type 06	100378
	GSII	Ultra Wide	Type 06	100378
	GSII	Wide	Type 06	100378
	GSIII	Mini	Type 05	100376
	GSIII	Regular	Type 06	100378
	SSI	Regular	Type 04	100350
	SSI	Wide	Type 04	100350
	SSII	Regular	Type 04	100350
	SSII	Wide	Type 04	100350
	SSIII	Regular	Type 04	100350
	SSIII	Wide	Type 04	100350
	TSIII CA	3.5	Type 05	100376
	TSIII SA	3.0	Type 11	100372
	TSIII SA	4.0	Type 47	100478
	TSIII SA	4.5	Type 47	100478
	TSIII SA	5.0	Type 47	100478
	TSIII SA Fixture – NoMount	4.5	Type 47	100478
	TSIII SA Ultra-Wide	6.0	Type 04	100350
	TSIII SA Ultra-Wide	7.0	Type 04	100350
	USII	Mini	Type 02	100355
	USII	Regular	Type 01	100353
	USII	Ultra Wide	Type 03	100357
	USII	Wide	Type 03	100357
	USIII	Mini	Type 02	100355
	USIII	Regular	Type 01	100353
	USIII	Ultra Wide	Type 03	100357
	USIII	Wide	Type 03	100357
	USIV Conical	Mini	Type 08	100366
<b>OSTEOPLANT</b>	Two-piece		Type 10	100370
<b>OSTEO-TI</b>	Combi Dental Implant	3.25	Type 26	100425
	Combi Dental Implant	3.75	Type 26	100425
	Combi Dental Implant	4.5	Type 26	100425
	Combi Dental Implant	5.5	Type 26	100425
<b>OTICON MEDICAL</b>	Ponto Implant	3.0	Type 09	100368
	Ponto Implant	4.0	Type 09	100368
	Ponto Implant	Wide	Type 09	100368

Implant Company	Model	Size	SmartPeg	Art. No
<b>PALTOP</b>	PALTOP	NP	Type 38	100455
	PALTOP	SP	Type 49	100480
<b>PHIBO</b>	Aurea	NP	Type 63	100554
	Aurea	RP	Type 27	100431
	Aurea	WP	Type 64	100555
	BNT	S3	Type 01	100353
	BNT	S4	Type 01	100353
	BNT	S5	Type 01	100353
	TSA	S3	Type 38	100455
	TSA	S4	Type 27	100431
	TSA	S5	Type 27	100431
	TSA Advance	S3	Type 38	100455
	TSA Advance	S4	Type 27	100431
	TSA Advance	S5	Type 27	100431
	TSH	S3	Type 01	100353
	TSH	S4	Type 01	100353
	TSH	S5	Type 01	100353
<b>P-I BRÄNEMARK PHILOSOPHY</b>	Internal Grip Amplified	3.5	Type 63	100554
	Internal Grip Amplified	4.3	Type 51	100482
	Internal Grip Amplified	5.1	Type 51	100482
	Internal Grip External Hex	3.5	Type 38	100455
	Internal Grip External Hex	3.75	Type 18	100398
	Internal Grip External Hex	4.1	Type 18	100398
	Morse Taper	3.3	Type 21	100402
	Morse Taper	3.75	Type 21	100402
Morse Taper	4.8	Type 21	100402	
<b>PLAN 1 HEALTH</b>	GTB	3.3	Type 57	100497
	GTB	3.6	Type 57	100497
	GTB	4.3	Type 57	100497
	Hexcel EHC, External Hex		Type 01	100353
	Hexcel IHC, Internal Hex		Type 26	100425
Hexcel IHT, Internal Hex		Type 26	100425	
<b>PLATON</b>	EL Implant Narrow	3.7	Type 60	100543
	EL Implant Regular	5.1	Type 47	100478
	Type I		Type 04	100350
	Type II		Type 04	100350
	Type IV Pro		Type 47	100478
<b>PRODENT ITALIA</b>	Perform SM	3.6	Type 49	100480
	Perform SM	4.0	Type 49	100480
	Perform SM	4.5	Type 07	100380
	Perform SM	5.0	Type 07	100380



Implant Company	Model	Size	SmartPeg	Art. No
<b>REVIOS</b>		3.8	Type 26	100425
		4.3	Type 26	100425
		5.0	Type 26	100425
<b>SDS</b>	SDS 2.0	3.8	Type 69	100589
<b>SERVO-DENTAL</b>	Trias	3.8	Type 06	100378
<b>SHINHUNG</b>	Luna	NP 3.5	Type 05	100376
	Luna	RP 4.0	Type 47	100478
	Luna	RP 4.5	Type 47	100478
	Luna	RP 5.0	Type 47	100478
	Sola	NP 3.5	Type 02	100355
	Sola	RP 4.0	Type 01	100353
	Sola	RP 4.5	Type 01	100353
	Sola	WP 5.0	Type 03	100357
	Stella	RP 4.0	Type 61	100544
	Stella	RP 4.5	Type 61	100544
	Stella	WP 4.5	Type 61	100544
	Stella	WP 5.0	Type 61	100544
<b>SIC IMPLANT SYSTEM</b>	Ace	3.4	Type 44	100466
	Ace	4.0	Type 44	100466
	Ace	5.0	Type 44	100466
<b>S.I.N IMPLANT SYSTEM</b>	External Hex	3.4	Type 01	100353
	External Hex	4.1	Type 01	100353
	External Hex	5.0	Type 01	100353
	Internal Hex	3.8	Type 15	100386
	Internal Hex	4.5	Type 15	100386
	Internal Hex	5.5	Type 15	100386
	Strong SW HI	4.5	Type 48	100479
	Strong SW Morse	4.5	Type 53	100484
<b>SOUTHERN IMPLANTS</b>	External Hex		Type 01	100353
	Extra-oral	3.75	Type 01	100353
	IT Connection		Type 04	100350
	MAX	7.0	Type 01	100353
	MAX	8.0	Type 01	100353
	MAX	9.0	Type 01	100353
	MAX-IT	7.0	Type 07	100380
	MAX-IT	8.0	Type 07	100380
	MAX-IT	9.0	Type 07	100380
	Tri-MAX	7.0	Type 24	100412
	Tri-MAX	8.0	Type 40	100462
Tri-MAX	9.0	Type 40	100462	
Tri-Nex	3.5	Type 12	100359	

Implant Company	Model	Size	SmartPeg	Art. No
<b>SOUTHERN IMPLANTS</b> (Continued)	Tri-Nex	4.3	Type 13	100361
	Tri-Nex	5.0	Type 13	100361
<b>STERNGOLD DENTAL</b>	External Hexagon	RP 3.75	Type 01	100353
	External Hexagon	RP 4.0	Type 01	100353
	External Hexagon	RP 5.0	Type 01	100353
	Stern IC Internal Octagon	RN 4.1	Type 04	100350
<b>STRAUMANN</b>	Bone Level Implants Roxolid SLA Active NC	3.3	Type 53	100484
	Bone Level Implants Roxolid SLA Active RC	4.1	Type 54	100485
	Bone Level Implants Roxolid SLA Active RC	4.8	Type 54	100485
	Bone Level Tapered Implants Roxolid SLA Active NC	3.3	Type 53	100484
	Bone Level Tapered Implants Roxolid SLA Active RC	4.1	Type 54	100485
	Bone Level Tapered Implants Roxolid SLA Active RC	4.8	Type 54	100485
	Standard Implants Roxolid SLA Active RN	3.3	Type 04	100350
	Standard Implants Roxolid SLA Active RN	4.1	Type 04	100350
	Standard Implants Roxolid SLA Active RN	4.8	Type 04	100350
	Standard Implants Roxolid SLA Active WN	4.8	Type 04	100350
	Standard Plus Implants Roxolid SLA Active NNC	3.3	Type 68	100560
	Standard Plus Implants Roxolid SLA Active RN	3.3	Type 04	100350
	Standard Plus Implants Roxolid SLA Active RN	4.1	Type 04	100350
	Standard Plus Implants Roxolid SLA Active RN	4.8	Type 04	100350
	Standard Plus Implants Roxolid SLA Active WN	4.8	Type 04	100350
	Tapered Effect Implants Roxolid SL Active RN	3.3	Type 04	100350
	Tapered Effect Implants Roxolid SL Active RN	4.1	Type 04	100350
	Tapered Effect Implants Roxolid SL Active WN	4.8	Type 04	100350
	Bone Level Implants Titanium SLA NC	3.3	Type 53	100484
	Bone Level Implants Titanium SLA RC	4.1	Type 54	100485
	Bone Level Implants Titanium SLA RC	4.8	Type 54	100485
	Standard Implants Titanium SLA RN	3.3	Type 04	100350
	Standard Implants Titanium SLA RN	4.1	Type 04	100350
	Standard Implants Titanium SLA RN	4.8	Type 04	100350
	Standard Implants Titanium SLA WN	4.8	Type 04	100350
	Standard Plus Implants Titanium SLA RN	3.3	Type 04	100350
	Standard Plus Implants Titanium SLA RN	4.1	Type 04	100350
	Standard Plus Implants Titanium SLA RN	4.8	Type 04	100350
	Standard Plus Implants Titanium SLA WN	4.8	Type 04	100350
	Tapered Effect Implants Titanium SLA RN	3.3	Type 04	100350
	Tapered Effect Implants Titanium SLA RN	4.1	Type 04	100350
	Tapered Effect Implants Titanium SLA WN	4.8	Type 04	100350
	Extra Oral		Type 01	100353
	Narrow Neck		Type 17	100390
	Ortho		Type 09	100368
	Palatal		Type 09	100368
<b>SWEDEN &amp; MARTINA</b>	EZ		Type 07	100380
	Global	3.8	Type 48	100479

Implant Company	Model	Size	SmartPeg	Art. No
<b>SWEDEN &amp; MARTINA</b> (Continued)	Global	4.3	Type 48	100479
	Global	4.8	Type 48	100479
	Global	5.5	Type 48	100479
	Kohno - new platform	3.8	Type 32	100440
	Kohno - old platform	3.8	Type 21	100402
	Kohno - new platform	4.25	Type 07	100380
	Kohno - old platform	4.25	Type 07	100380
	Kohno - new platform	5.0	Type 07	100380
	Kohno - old platform	5.0	Type 07	100380
	Premium - new platform	3.3	Type 32	100440
	Premium - old platform	3.3	Type 21	100402
	Premium - new platform	3.8	Type 32	100440
	Premium - old platform	3.8	Type 21	100402
	Premium- new platform	4.25	Type 07	100380
	Premium- old platform	4.25	Type 07	100380
	Premium -new platform	5.0	Type 07	100380
	Premium - old platform	5.0	Type 07	100380
<b>SYBRON</b>	Endopore External	3.5	Type 30	100436
	Endopore External	4.1	Type 01	100353
	Endopore External	5.0	Type 01	100353
	Entegra External	3.25	Type 01	100353
	Entegra External	4.0	Type 01	100353
	Pitt-Easy	4.0	Type 24	100412
	Pitt-Easy	4.9	Type 24	100412
	SybronPRO XRT	3.3	Type 06	100378
	SybronPRO XRT	3.5	Type 06	100378
SybronPRO XRT	4.1	Type 06	100378	
<b>TA-DENT</b>		3.3	Type 38	100455
		4.8	Type 06	100378
		5.6	Type 49	100480
<b>TBR</b>	External Hexagon		Type 01	100353
	Internal Octagon		Type 18	100398
	M-implant	3.2	Type 21	100402
	M-implant	3.9	Type 21	100402
	M-implant	4.7	Type 21	100402
	Z1-M implant	3.2	Type 21	100402
	Z1-M implant	3.9	Type 21	100402
Z1-M implant	4.7	Type 21	100402	
<b>TEKKA</b>	Inkone	3.5	Type 49	100480
	Inkone	4.5	Type 49	100480
	Krestal	3.5	Type 06	100378
	Krestal	4.0	Type 06	100378
	Krestal	4.2	Type 06	100378

Implant Company	Model	Size	SmartPeg	Art. No
<b>TEKKA (Continued)</b>	Krestal	5.1	Type 06	100378
	Progress	3.4	Type 06	100378
	Progress	4.0	Type 06	100378
	T-quest	4.0	Type 49	100480
	T-quest	4.5	Type 49	100480
	Twinkone	3.5	Type 17	100390
<b>TFI SYSTEM</b>	TFI	3.3	Type 26	100425
	TFI	3.75	Type 26	100425
	TFI	4.2	Type 26	100425
	TFI	5.0	Type 26	100425
<b>THOMMEN MEDICAL</b>	SPI System Platform	3.5	Type 22	100404
	SPI System Platform	4.0	Type 21	100402
	SPI System Platform	4.5	Type 15	100386
	SPI System Platform	5.0	Type 15	100386
	SPI System Platform	6.0	Type 15	100386
<b>TIDAL SPIRAL</b>		3.75	Type 27	100431
		4.2	Type 27	100431
		5.0	Type 27	100431
		6.0	Type 27	100431
<b>TRI DENTAL IMPLANTS</b>	TRI-Octa 4.8	3.75	Type 61	100544
	TRI-Octa 4.8	4.1	Type 61	100544
	TRI-Octa 4.8	4.7	Type 61	100544
	TRI-Vent	3.5 platform	Type 32	100440
<b>U-IMPL</b>	External Hex		Type 01	100353
	Internal Hex		Type 32	100440
<b>WARANTEC</b>	Oneplant Implant	4.3	Type 24	100412
	Oneplant Internal Tapered		Type 06	100378
<b>WIELAND IMPLANT SYSTEM</b>		3.5	Type 45	100472
		4.3	Type 45	100472
		5.0	Type 45	100472
<b>ZIMMER DENTAL</b>	AdVent	4.5	Type 35	100446
	AdVent	5.7	Type 34	100444
	Screw-vent 3.5	3.3	Type 32	100440
	Screw-vent 3.5	3.7	Type 32	100440
	Screw-vent 4.5	4.7	Type 27	100431
	Tapered Screw-vent 3.5	3.7	Type 32	100440
	Tapered Screw-vent 3.5	4.1	Type 32	100440
	Tapered Screw-vent 4.5	4.7	Type 27	100431

Implant Company	Model	Size	SmartPeg	Art. No	
<b>ZIMMER DENTAL</b> (Continued)	Tapered Screw-vent	5.7	6.0	Type 33	100442
	Spline	3.25	3.25	Type 30	100436
	Spline	3.75	3.75	Type 30	100436
	Spline	5.0	5.0	Type 30	100436
	SwissPlus	3.8	3.8	Type 31	100438
	SwissPlus	4.8	4.8	Type 04	100350

<b>ZUGA MEDICAL</b>	Bone Level	5.0	5.0	Type 65	100557
---------------------	------------	-----	-----	---------	--------

Abutment System	Model	Size	SmartPeg	Art. No	
<b>COCHLEAR</b>	Baha BA 210	5.5	5.5	Type 55	100487
	Baha BA 210	8.5	8.5	Type 55	100487
	Baha BA 300, BIA 300	6.0	6.0	Type 55	100487
	Baha BA 300, BIA 300	9.0	9.0	Type 55	100487
	Baha BA 400, BIA 400	6.0	6.0	Type 55	100487
	Baha BA 400, BIA 400	8.0	8.0	Type 55	100487
	Baha BA 400, BIA 400	10.0	10.0	Type 55	100487
	Baha BA 400, BIA 400	12.0	12.0	Type 55	100487
	Baha BA 400, BIA 400	14.0	14.0	Type 55	100487
	Vistafix Standard abutment	3.0	3.0	Type A3	100363
	Vistafix Standard abutment	4.0	4.0	Type A3	100363
	Vistafix Standard abutment	5.5	5.5	Type A3	100363
	Vistafix Standard abutment	7.0	7.0	Type A3	100363
	Vistafix VXA 300	3.5	3.5	Type A3	100363
	Vistafix VXA 300	4.5	4.5	Type A3	100363
	Vistafix VXA 300	6.0	6.0	Type A3	100363
	Vistafix VXA 300	7.5	7.5	Type A3	100363

<b>NEOSS</b>	Access abutment			Type 52	100483
--------------	-----------------	--	--	---------	--------

<b>NOBEL BIOCARE</b>	Brånemark system MUA	NP		Type A3	100363
	Brånemark system MUA	RP		Type A3	100363

<b>OTICON MEDICAL</b>	Ponto Implant abutment	4.0		Type 55	100487
-----------------------	------------------------	-----	--	---------	--------

<b>THOMMEN MEDICAL</b>	SPI System Bar abutment			Type 36	100447
	VARIOMulti abutment			Type 25	100415

# SmartPeg

## **A fresh start every time!**

To avoid damage to the implant, the SmartPeg material is quite soft, which limits the lifetime of the threads. Therefore, SmartPegs are single use items.





